

"Investing in Africa's Future"

COLLEGE OF HEALTH, AGRICULTURE & NATURAL SCIENCES DEPARTMENT OF BIOMEDICAL AND LABORATORY SCIENCES

NSLS203 Haematology Practical

END OF SECOND SEMESTER EXAMINATIONS

April/May 2023

LECTURER: Mrs E. Govore

DURATION: 3 HOURS

INSTRUCTIONS

- 1. Write your candidate number on the space provided on top of each page
- 2. Answer all questions in sections A, B and C
- 3. Credit will be given for logical, systematic and neat presentations

Candidate number.....

Ouestion 1

- a) Using the slide making and staining procedure given below make a thin smear of the provided sample **G** and stain it with the given Romanowsky stain (10 marks)
- **b)** List 5 types of Romanowsky stains. (5)
- c) What are the different parts of the thin smear (6)
- d) Which is the best area to conduct a differential count (2)
- e) There are 5 types of blood smears. List the different types of blood smears (2)
- f) Viewing the smear, give general comments on red cells, platelets and white cells (10)

Procedure for slide making

Thin smear

STEP ONE:

- Place clean glass slide on a flat surface
- · Add one small drop of blood to one end

STEP TWO:

- · Take another clean slide
- holding at an angle of about 45 deg
- touch the blood with one end of the slide so the blood runs along the edge of the slide by capillary action
- Push carefully along the length of the first slide to produce a thin smear of blood

STEP THREE:

• Make 2 smears, allow to air dry, and label clearly with your candidate number

STEP FOUR:

•Fix the thin film with methanol (100% or absolute) and let it dry completely before staining

Procedure for thin film staining

- 1. Fresh working Giemsa stain has been prepared for you
- 2. Place slides onto the staining rack and add the working Giemsa stain
- 3. Stain for 10minutes
- 4. Rinse slides with Giemsa buffer
- 5. Dry the slides upright in a rack.

(andidate)	numhar		
Calluluate	HUHHDEL	 	

Question 2

You are provided with sample ${\bf C}$ and you are required to perform the test indicated in the procedure below

Procedure for the question

- **1.** Prewarm reagents and patient sample specimens to 37°C prior to performing this test.
- **2.** Pippet 100ul of Citrated plasma into a labeled 12X75mm test tube.
- 3. Add 100ul of R1 to the test tube.
- **4.** Mix and incubate for five minutes at 37°C.
- **5.** Add 100ul of R2 to the test tube.
- **6.** On addition of R2 start stop watch simultaneously.
- **7.** Mix and determine coagulation time by tilting tube back and forth (checking for clot formation)
- **8.** Stop stopwatch as soon as clot forms
- 9. Record and report results.

Questions:

- (a) Give the name of the test which you have performed and explain the principle of this test. Which coagulation pathway is being measured in your test? (10)
- (b) State and discuss the results obtained using this method (15)
- (c) Given that you conducted Prothrombin (PT) and Activated Partial Thromboplastin time(APTT) and the following results were obtained state whether intrinsic, extrinsic or common pathway is affected. P = prolonged, N= Normal (5)

test	PT	APTT	Pathway affected
i.	N	P	
ii.	P	P	
iii.	P	N	
iv.	N	N	

Question 3 (30 MARKS)

1. State the different cells or images as projected on the microscope or picture from number 1-20