NSLS104/JANJUN22



# COLLEGE OF HEALTH, AGRICULTURE AND NATURAL SCIENCES

# DEPARTMENT OF BIOMEDICAL AND LABORATORY SCIENCES

#### **NSLS104: CLINICAL PATHOLOGY PRACTICAL**

#### END OF SECOND SEMESTER FINAL EXAMINATIONS

#### APRIL/MAY 2022

# LECTURER: MR G. MALUNGA

## 22 APRIL 2022 0900 HRS

# **DURATION: 3 HOURS**

# **INSTRUCTIONS**

Answer **all** questions on the separate answer sheet provided.

Mark allocation for each question is indicated at the end of the question

Credit will be given for logical, systematic and neat presentations

1.0 ml

\_

1.0 ml

1.0 ml

#### **Question** 1

A woman who is about 8 months pregnant visits an Antenatal Clinic complaining of lower abdominal pain and swollen lower limbs. The doctor asked her to produce a urine sample for protein measurement. You are required to determine the amount of protein in the urine sample labeled M from the woman using the following procedure.

#### Urine protein determination procedure

- 1. Using pH papers, test the pH of the urine. If it is neutral or alkaline, add 3 drops of glacial acetic acid to it.
- 2. Prepare a pair of test tubes for the standard and a pair for the sample **M** as shown in Table 1 below. Label one tube B (for blank) and the other one T (for the test). The provided protein standard solution has a protein concentration of 35 g/l

| iable 1. i reparation of reaction tubes |                            |                  |        |                       |        |
|---|----------------------------|------------------|--------|-----------------------|--------|
|   |                            | PROTEIN STANDARD |        | URINE SAMPLE <b>M</b> |        |
|   |                            | BLANK            | TEST   | BLANK                 | TEST   |
|   | SULPHOSALICYLIC ACID (SSA) | -                | 1.5 ml | -                     | 1.5 ml |
|   | DEIONISED WATER            | 1.5 ml           | -      | 1.5 ml                | -      |

#### Table 1: Prenaration of reaction tubes

- 3. Mix the contents of each tube thoroughly and leave them for 5 minutes.
- 4. Remix the pair of the protein standard tubes and immediately set the spectrophotometer at zero absorbance with the B tube at 590 nm. Read and record the absorbance of the corresponding T tube.

1.0 ml

- 5. Repeat step 4 for the pair of the urine sample tubes.
- 6. Carry out a semi-quantitative determination of the protein in urine **M** using the provided urine dipsticks and record your results.

#### Questions

URINE SAMPLE **M** 

PROTEIN STANDARD

- a) Why is alkaline urine not suitable for the SSA method? [2]
- b) Record your absorbance results for the protein standard and urine sample in the form of a table. [20]
- c) Calculate the concentration of protein in urine  $\mathbf{M}$ . [5]
- d)Record your protein results from the semi-quantitative determination of the urine protein. [5]
- e) Explain your results in relation to the symptoms presented by the woman. [3]
- f) Which other laboratory tests can be done on the urine in order to come up with a definite diagnosis? [5]

# **TOTAL MARKS: 40**

## Question 2

You are required to carry out a microscopy examination and biochemical analysis of a urine sample labeled  $\mathbf{U}$  from a 55 year-old woman complaining of painful urination. Use the following procedure.

## Procedure

- 1. Pour the provided urine into a conical centrifuge tube to about two thirds full.
- 2. Note the appearance of the urine.
- 3. Dip a urine reagent strip into the urine and record the biochemical findings.
- 4. Centrifuge the urine in a centrifuge at 2500 rpm for 5 minutes.
- 5. Decant the supernatant completely.
- 6. Resuspend the sediment by tapping the bottom of the tube and transfer one drop of the sediment onto a microscope glass slide and cover with a cover slip.
- 7.Examine the wet preparation using a microscope starting with the 10X objective and then move on to the 40X objective.
- 8. Record your findings.

## Questions

- a) Record the appearance of the urine. [1]
- b) Record the following biochemical findings of the urine. [22]

Leucocytes Nitrite Urobilinogen Protein pH Blood Specific Gravity Ketones Bilirubin

Glucose

c) Record the following microscopy findings of the urine. [25]

WBC RBC Epithelial cells Yeasts Casts Crystals S. haematobium Others

d) What diagnosis can you make from your results? Support your diagnosis with your results? [12]

# **TOTAL: 60 MARKS**