

"Investing in Africa's Future"

COLLEGE OF HEALTH, AGRICULTURE AND NATURAL SCIENCES DEPARTMENT OF HEALTH SCIENCES BACHELOR OF MEDICAL LABORATORY SCIENCES HONOURS

NSLS101: LABORATORY PRINCIPLES PRACTICAL

END OF FIRST SEMESTER EXAMINATIONS

JUNE 2019

LECTURER: MR G. MALUNGA

DURATION: 3 HOURS

INSTRUCTIONS

Answer all questions in both sections A and B on separate answer sheets provided.

Mark allocation for each question is indicated at the end of the question.

Credit will be given for logical, systematic and neat presentations.

SECTION A: SPOT EXAM [40MARKS]

There are 10 stations which have been set up in the laboratory numbered 1 to 10. For each station (a) Name the item/s shown

(b) State it's/their use in the laboratory

SECTION B [60 MARKS]

Question 1

You are provided with Sodium Chloride granules, 3.5% Sodium Hypochlorite solution and absolute ethanol. Use these provided chemicals to prepare the following solutions. Fully describe the method of preparation of each solution showing all the calculations, masses and volumes used. Label your solutions clearly showing your candidate number and the name of the solution. Submit the prepared solutions for marking.

- (a) 75ml of physiologic saline.[10]
- **(b)** 90ml of 1% Sodium hypochlorite solution using 3.5% Sodium hypochlorite solution.[10]
- **(c)** 80ml of 70% ethanol.[10]
- **(d)** State one laboratory use of each of the solutions which you have prepared. [3]
- (e) Describe how a balance can be maintained in good working condition. [2]

TOTAL: 35 MARKS

Question 2

You are provided with a urine sample labeled **W.** Mix the urine thoroughly and pour it in a centrifuge tube and centrifuge it at 3000 rpm for 5 minutes. Pour away the supernatant and keep the sediment. Tap off the sediment and place a sufficient amount of the sediment on a glass slide and cover with a cover slip. Observe the sediment under a microscope.

- (a) Obtain a clear focus of one field and make a good drawing of what you are observing under the X40 objective lens. [20]
- (b) Which safety precaution measures must be considered when using a
 - (i) Microscope
 - (ii) Centrifuge [5]

TOTAL: 25 MARKS