

"Investing in Africa's future" COLLEGE OF HEALTH, AGRICULTURE AND NATURAL SCIENCES

DEPARTMENT OF BIOMEDICAL AND MEDICAL LABORATORY SCIENCES

BACHELOR OF MEDICAL LABORATORY SCIENCES HONOURS

NSLS211: HISTOLOGY PRACTICAL

END OF SEMESTER FINAL EXAMINATIONS

NOVEMBER 2023

LECTURER: DR M. SALISSOU

DURATION: 3 HOURS

INSTRUCTIONS

SECTION A: Spot exam: Answer all questions on separate sheets provided

SECTION B: Answer all questions on separate answer sheets provided

SECTION C: Answer all questions on separate answer sheets provided

SECTION A

Spot exam: Answer all questions.

- 1. A 1, A2, A3, A4 and A5 are special stains in histology. For each of the stains
 - a) Identify the stain
 - b) report the results
 - c) justify your answers

[50]

SECTION B

Clinical case: Answer all questions

1. A 45-year-old man presents with increasing abdominal girth and yellow discoloration of his skin and sclera. Physical examination reveals hepatomegaly and jaundice.. a Histological laboratory procedure for the identification was performed and a stain of a liver biopsy is shown in Fig 1. You have been provided with the prepared slide *TB1* for examination. What is the major intracellular iron storage protein in this patient's hepatocytes?




Fig 1

- a) Identify the type of stain:
- b) Carefully examine the slide and report your findings (5).
- c) Explain the principle of the test (10)

d) If your tissue was a liver biopsy explain all process to obtain paraffin embedded block of the tissue (10)

SECTION C

Answer all questions

1. Carry out the following staining protocol on the section provided using the method below (25).

Method

i. Dewax sections and bring sections to water as follow:

a.	Xylene	10 dips
b.	Xylene	10 dips
c.	Xylene	10 dips
d.	Absolute alcohol	10 dips
e.	Absolute alcohol	10 dips
f.	95% alcohol	10 dips
g.	85% alcohol	10 dips
h.	75% alcohol	10 dips
i.	Rinse in tap water	1 minute

- j. Place slide in ammoniacal silver solution and incubate at 58 to 60 degree in water bath for 30 to 60 minutes until section become yellowish brown.
- k. Rinse in distilled water
- 1. Place the slide in1% neutral red for two minutes
- m. Wash in distilled water
- n. Dehydrate and clear the section as follows:
- ii. 75% alcohol 10 dips
- iii. 85 % alcohol 10 dips
- iv. 95% alcohol 10 dips
- v. Absolute alcohol 10 dips
- vi. Absolute alcohol 10 dips
- vii. Xylene 10 dips.
- viii. Xylene 10 dips
- ix. Mount in poly-x mountant
- x. Label your section using provided stickers.
- a) What is the name of this stain? (5)
- b) What is the principle of this stain? (15)
- c) What further precautions to ensure (not mentioned above) can help identify substances in the tissue (5)